

FloorBed™ 1






FloorBed™ 1


TECHNICAL SPECIFICATION

Dimensions:	36" x 89"	Lifting range:	2.75" - 25" (excl. mattress)
Platform size:	36" x 79"	Profiling functions:	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Safe working load:	408 lbs		
Max patient weight:	330 lbs	Part number:	NSB-0-FL1-200US




Why FloorBed™?

NURSING AND LOW BEDS



FLOORBED™



FloorBed™ 1 Features



Width adjustable mattress guides



No trailing wires



Full nursing height



Concealed wheels



Backrest and kneebreak functions



Mattress platform height of just 2.75 inches – similar to a tennis ball!

FloorBed™ 1 Accessories



Assist bar to assist with bed mobility



Trapeze to assist with bed mobility